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
图 4. Figure 4. Resolution vs. Input Counts Rate (ICR) for Various Peaking Times for FAST SDD.

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
General	
Detector Type	Silicon Drift Detector (SDD) with CMOS preamplifier
Detector Size	25 mm ²
Silicon Thickness	500 μm, See efficiency curves
Collimator	Internal MultiLayer Collimator (ML)
Energy Resolution @ 5.9 keV (55Fe)	125 - 140 eV FWHM depending on peaking time
Peak to Background	20,000:1 (ratio of counts from 5.9 keV to 1 keV) (typical)
Detector Window Options	Beryllium (Be): 0.5 mil (12.5 μm) or 0.3 mil (8 μm)
	C Series Low energy windows
Charge Sensitive Preamplifier	Amptek custom reset preamplifier
Gain Stability	<20 ppm/°C (typical)
XR-100FastSDD Case Size	3.00 x 1.75 x 1.13 in (7.6 x 4.4 x 2.9 cm)
XR-100FastSDD Weight	4.4 ounces (125 g)
Total Power	<2 Watt
Warranty Period	1 Year
Typical Device Lifetime	5 to 10 years, depending on use
Operation conditions	0°C to +50°C
Storage and Shipping	Long term storage: 10+ years in dry environment
	Typical Storage and Shipping: -20°C to +50°C, 10 to 90% humidity non condensing
	TUV Certification
	Certificate #: CU 72072412 02
	Tested to: UL 61010-1: 2004 R7 .05
	CAN/CSA-C22.2 61010-1: 2004
Inputs	
Preamplifier Power	±8 to 9 V @ 15 mA with no more than 50 mV peak-to-peak noise
Detector Power	-70 to -200 V @ 25 μA very stable <0.1% variation
Cooler Power	Current = 450 mA maximum, voltage = 3.5 V maximum with <100 mV peak-to-peak noise

	Note: the XR-100SDD includes its own temperature controller
Outputs	
Preamplifier Sensitivity	4 mV/keV typical (may vary for different detectors)
Preamplifier Polarity	Positive signal output (1 kohm maximum load)
Preamplifier Feedback	Reset
Temperature Monitor Sensitivity	PX5: direct reading in Kelvin through software
Options	
X-123FastSDD	 The silicon drift detector (FastSDD) is also available in the X-123SDD Silicon Drift Detector System configuration. The X-123SDD configuration includes the detector, preamplifier, DP5 digital pulse processor and MCA, and the PC5 power supply. All that is needed is a +5 Volts DC input and a USB, RS232, or Ethernet connection to your computer. X-123SDD X-Ray Detector System
Vacuum Accessories	The FastSDD is compatible with all Amptek vacuum accessories
OEM	The FastSDD is compatible with all Amptek OEM configurations

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