


# HeNe Laser with polarization maintaining fiber optics

## Universal Beam Source for Interferometry and Frequency Standards

High coupling efficiency, transport stability, and polarization maintenance are characteristics of **Schäfter + Kirchhoff** helium-neon (HeNe) lasers with fiber optics. HeNe lasers with wavelength 543.5, 594.1, and 632.8 nm are available with output

powers at the fiber ends ranging from approx. 1 mW (all wavelengths) to 30 mW (632.8 nm only). For low-maintenance laser operation in rough environments, **Schäfter + Kirchhoff** provides shock and vibration absorbing mounting consoles.



**Some fiber coupled HeNe lasers available from Schäfter + Kirchhoff:**

- A** HeNe laser 632.8 nm, output power ex fiber 30 mW, with mechanical attenuator, see page 4
- B** Frequency stabilized HeNe laser 632.8 nm, with Faraday isolator, mounted on shock and vibration absorbers, see page 5.
- C** HeNe laser 632.8 nm, 0.7 mW ex fiber, in mounting console with electro-magnetic shutter (option)
- D** HeNe laser 632.8 nm, 13 mW ex fiber, see page 3.
- E** Same laser as above, with mechanical attenuator and shock and vibration absorbing mounting console.


See following pages for selection of laser source and mounting option.

**Laser Beam Analysis**  
24.01.2007  
Schäfter & Kirchhoff Hamburg

HeNe laser type: HMP 515  
Wavelength:  $\lambda = 632.8$  nm  
Fiber: PMC-638-4,1-100-3-AP-158-P  
Cut-off:  $\lambda_{co} = 542$  nm  
MFD = 4.1  $\mu$ m  
NR = 0.12

| Wavelength (nm)            | 543.5 | 594.1 | 632.8 | 632.8 |
|----------------------------|-------|-------|-------|-------|
| Output Power ex Fiber (mW) | 0.7   | 1.4   | 0.7   | 1-30  |
| Lin. polarized             | yes   | yes   | yes   | yes   |
| Freq. stabilized           | -     | -     | yes   | -     |

### Applications



## HeNe Laser with singlemode fiber

### All-round beam source for interferometry and frequency standards

For more than 40 years HeNe lasers are the base of interferometers, spectroscopic, and holographic measurement setups. If wavelength and frequency stability are required, HeNe lasers are superior over many other laser sources with respect to costs, live-time, and reliability.

For HeNe lasers with polarization maintaining fiber optics, there is an increased number of possible applications in research and industry. Coupling the laser beam to a singlemode fiber cable makes a spatial separation of laser beam source and measurement setup possible.

Voluminous, heat dissipating lasers can be placed in a safe distance from the experiment. An additional advantage there is from a laser safety point of view. Even in large setups, it is possible to confine all laser beams to optical fibers, and to completely avoid free space beams.

**Made in Germany**

## Contents

### HeNe Lasers

#### HeNe lasers red

Page 3


- Wavelength 632.8 nm
- Fiber-coupled to singlemode or polarization maintaining fiber cable
- Output power 1.5 - 12 mW ex fiber
- Mounting console with shock absorbers
- Optionally mechanical attenuator



#### HeNe lasers green, yellow, and infrared

Page 3

- Wavelengths 543.5 nm (green), 594.1 nm (yellow), and 1523 nm (infrared)
- Fiber-coupled to singlemode or polarization maintaining fiber cable
- Output power 0.7 mW (green), 1 mW (yellow), and 0.6 mW (infrared) ex fiber
- Mounting console with shock absorbers
- Optionally mechanical attenuator



#### High power HeNe laser 633 nm - output power 28 mW ex fiber

Page 4

- Wavelength 632.8 nm
- Fiber-coupled to polarization maintaining fiber cable
- Output power 28 mW ex fiber
- Optionally mechanical attenuator



#### Compact HeNe laser 633 nm / 0.6 mW ex fiber - standard and industrial version

Page 4

- Fiber-coupled to singlemode or polarization maintaining fiber cable
- Laser class 2 for adjustment purposes
- Mounting console with shock absorbers
- Optionally mechanical or electrical shutter or mechanical attenuator



#### Frequency stabilized HeNe laser 633 nm

Page 5

- Wavelength 632.8 nm
- Output power 0.7 mW ex fiber
- frequency stabilized (typ. < 500 kHz over some minutes, < 5 MHz over some hours)
- Fiber coupled to polarization maintaining fiber cable
- Faraday isolator avoiding back reflections into the laser resonator
- Mounting console with shock absorbers
- Optionally mechanical attenuator




### Components and Accessories

#### Components for do-it-yourself laser coupling

Page 6


- Laser beam couplers 60SMS-...
- Adapters and attenuators 60A19.5-F-...
- Fiber cables polarization maintaining PMC-... or singlemode SMC-...
- Fiber collimators 60FC-...



#### Mounting consoles and accessories

Page 6


- With shock and vibration absorbing elements
- Mounting console with integrated coupling flange for an increased coupling stability



#### Custom coupling of HeNe lasers

Page 6

- Fiber-coupling of customer supplied HeNe laser, example of coupling a 10 years old laser



#### Dimensional drawings for power supplies

Page 7

- All fiber-coupled HeNe lasers are shipped with power supply



#### Laser safety

Page 8

- Laser safety goggles
- IR detection cards




## Overview: HeNe lasers 633, 543, 594, and 1523 nm - output power 0.6 - 28 mW ex fiber

| Table 1 HeNe Laser |                 |                     |              |         |                |                    |            |                        |                   |           |                     |                   |                  |             |                    |      |                           |
|--------------------|-----------------|---------------------|--------------|---------|----------------|--------------------|------------|------------------------|-------------------|-----------|---------------------|-------------------|------------------|-------------|--------------------|------|---------------------------|
| Column             | 1               | 2                   | 3            | 4       | 5              | 6                  | 7          | 8                      | 9                 | 10        | 11                  | 12                | 13               | 14          | 15                 | 16   | 17                        |
| No                 | Wavelength [nm] | Power ex fiber [mW] | Polarization | Type    | Supply Voltage | Shutter/Attenuator | Fiber Type | Fiber Connector Option | Fiber Length [cm] | Shutter * | Laser Diameter [mm] | Laser Length [mm] | Drawing see Page | Laser Class | Power Supply Style | NA   | MFD 1/e <sup>2</sup> [μm] |
| Red                |                 |                     |              |         |                |                    |            |                        |                   |           |                     |                   |                  |             |                    |      |                           |
| 1                  | HeNe 633        | 0.6                 | P            | 1107 P  | ..             | ..                 | ..         | ..                     | ..                | -         | 31.6                | 178               | 4                | 2           | A                  | 0.12 | 4.1                       |
| 2                  | HeNe 633        | 1.5                 | P            | LHP 121 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 272               | 3                | 3B          | B                  | 0.12 | 4.1                       |
| 3                  | HeNe 633        | 3.5                 | P            | LHP 151 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 396               | 3                | 3B          | B                  | 0.12 | 4.1                       |
| 4                  | HeNe 633        | 7                   | P            | LHP 991 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 484               | 3                | 3B          | C                  | 0.12 | 4.1                       |
| 5                  | HeNe 633        | 12                  | P            | LHP 925 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 637               | 3                | 3B          | C                  | 0.12 | 4.1                       |
| 6                  | HeNe 633        | 28                  | P            | LHP 928 | ..             | ..                 | ..         | ..                     | ..                | x         | rect.               | 1030              | 4                | 3B          | C                  | 0.12 | 4.1                       |
| Green              |                 |                     |              |         |                |                    |            |                        |                   |           |                     |                   |                  |             |                    |      |                           |
| 7                  | HeNe 543        | 0.7                 | P            | LGP 121 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 510               | 3                | 3B          | B                  | 0.12 | 3.5                       |
| Yellow             |                 |                     |              |         |                |                    |            |                        |                   |           |                     |                   |                  |             |                    |      |                           |
| 8                  | HeNe 594        | 1.0                 | P            | LYP 173 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 456               | 3                | 3B          | B                  | 0.12 | 3.9                       |
| Infrared           |                 |                     |              |         |                |                    |            |                        |                   |           |                     |                   |                  |             |                    |      |                           |
| 9                  | HeNe 1523       | 0.6                 | P            | LIP 171 | ..             | ..                 | ..         | ..                     | ..                | x         | 44.5                | 456               | 3                | 3B          | C                  | 0.13 | 9.1                       |

|                                                                                                                                                                                                                                                                                                                                                 |                   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|
| HeNe - 633 - 3.5 - P - LHP151 - 2 - A - 28 - 0 - 150                                                                                                                                                                                                                                                                                            | <b>Order Code</b> |
| Wavelength [nm]                                                                                                                                                                                                                                                                                                                                 |                   |
| Optical power ex fiber [mW]                                                                                                                                                                                                                                                                                                                     |                   |
| Polarization:<br>linear (standard) ..... P<br>random (option) ..... R                                                                                                                                                                                                                                                                           |                   |
| Laser type                                                                                                                                                                                                                                                                                                                                      |                   |
| Supply voltage:<br>230V (standard) ..... 2<br>110 V ..... 1<br>12 V ** ..... 3                                                                                                                                                                                                                                                                  |                   |
| Shutter/Attenuator<br>mechanical shutter * ..... S<br>mechanical attenuator ..... A<br>w/o ..... 0                                                                                                                                                                                                                                              |                   |
| Fiber / connector type:<br>10 = singlemode fiber cable, FC-PC connector (0° polish)<br>18 = singlemode fiber cable, FC-APC connector (8° polish)<br>20 = PM singlemode fiber cable, FC-PC connector (0°-polish)<br>28 = PM singlemode fiber cable, FC-APC connector (8°-polish), standard<br>Connector type ST, DIN-Avio, and E-2000 on request |                   |
| Fiber length [cm]                                                                                                                                                                                                                                                                                                                               |                   |
| Connector option:<br>0 = standard<br>C = core centered < 0.25 μm (singlemode only)                                                                                                                                                                                                                                                              |                   |


\* In case the laser has its own shutter, see table, column 11, an external shutter is dispensable  
 \*\* 12 V OEM power supply only available for laser type 1107P, Table 1 row 1



**HeNe laser**, wavelength 632.8, 545.5, 543.5, 594.1, or 1523 nm (Tab. 1, No. 2 - 6 and 7 - 9).

with laser beam coupler and polarization maintaining singlemode fiber cable. A diameter of 44.5 mm is common to these lasers, the tube length varies with laser power and wavelength. The **laser beam coupler 60SMS-...** is held by the adapter **60A19,5-F-AT** (with integrated attenuator, as shown), or the adapter **60A19,5-F** (without attenuator, see p. 6).

For mounting of the adapter and the laser, several options are available (see p. 6).




- Simple configuration:** four screws are used to fasten the adapter to the laser's front plate
- Same as before, additionally, mounting consoles MC-MG-44,5-R** with elastomer shock mounts (set of two). Shock and vibrations are damped, thermal deformation of the laser housing, which would be present in case of a rigid mounting of the laser, is prevented.

- Optimum configuration for maximum stability:** a metallic mounting console MC-MG-44,5-F-S holds both laser and adapter. The mounting console rests on wire rope isolators, thrusts and accelerations are absorbed in all directions.




## High power HeNe laser 633 nm - output power 28 mW ex fiber



High power **HeNe 632.8 nm** laser with polarization maintaining singlemode fiber. Up to 30 mW laser power ex fiber.

The laser is shown with the optionally available **mechanical attenuator 60A19,5-F-AT**. Due to the singlemode fiber, adjustment of the laser power with the mechanical attenuator has no influence on the shape of the emerging Gaussian beam profile.



HeNe - 633 - 28 - P - LHP 928 - 2 - A - 28 - 0 - 300 **Order Code**

Supply voltage:

230V (standard) ..... 2

110V ..... 1

Attenuator ..... A

w/o ..... 0

Fiber length [cm]

Fiber / connector type:



20 = PM singlemode fiber cable, FC-PC connector

(0°-polish)

28 = PM singlemode fiber cable, FC-APC connector

(8°-polish), standard

## Compact HeNe laser 633 nm/ 0.6 mW ex fiber - standard and ruggedized industrial version



HeNe laser 632.8 with polarization maintaining singlemode fiber. Laser power ex fiber typically > 0.6 mW.

In the standard version (left) a special adapter provides a stable connection of laser beam coupler to the laser housing. On the right, a ruggedized industrial version with an integrated electro-magnetic shutter is shown, mounted on a vibration absorbing console.

\* Mechanical shutter, electrical shutter and attenuator available with mounting console only

HeNe - 633 - 0,7 - P - 1107 P - 2 - A - 28 - 0 - 150 **Order Code**

Polarization:

linear (standard) ..... P

random (option) ..... R

Supply voltage:

230V (standard) ..... 2

110 V ..... 1

230 V OEM ..... 3

110 V OEM ..... 4

12 V OEM ..... 5

Shutter/Attenuator

mechanical shutter \* ..... S

electrical shutter \* ..... E

mechanical attenuator \* ..... A

w/o ..... 0

Fiber length [cm]

Connector option:

0 = standard

C = core centered < 0.25 µm

(singlemode only)

Fiber / connector type:

10 = singlemode fiber cable, FC-PC connector (0° polish)

18 = singlemode fiber cable, FC-APC connector (8° polish)

20 = PM singlemode fiber cable, FC-PC connector

(0°-polish)

28 = PM singlemode fiber cable, FC-APC connector


(8°-polish), standard

Connector type ST, DIN-Avio, and E-2000 on request

## Frequency Stabilized HeNe Laser with polarization maintaining fiber optics

In frequency stabilized HeNe lasers a difference signal of two adjacent, orthogonally polarized modes is used for thermally correcting the resonator length and so compensating frequency variations. After selecting one of these modes by means of a polarizer, at output there is radiation with high polarization extinction and high frequency stability (typ. < 500 kHz over some minutes, < 5 MHz over some hours).

Schäfter+Kirchhoff couples frequency stabilized HeNe laser with high coupling efficiency (>75%) to polarization maintaining singlemod fiber cables. Since even marginal back scattering from a measurement object or from the laser beam coupler affects the control loop of the laser, a Faraday isolator (optical diode) is placed in between laser and fiber coupling.



- Wavelength 632.8 nm
- Laser power ex fiber >0.7 mW
- Polarization extinction at fiber end > 200:1
- Switchable between operation mode constant frequency and constant power
- Faraday isolator for protecting the control loop of disturbing back reflections, isolation > 30 dB, insertion loss > 0.5 dB

- Mechanical attenuator for reproducibly reducing the output power
- Polarization maintaining fiber type PANDA, mode field diameter 4.1 mm, numerical aperture 0.12
- Fiber cable Ø 3 mm with Keflar strain relief
- Fiber connector type FC-APC (8° polish)
- option: fiber cable with vacuum feed-through
- Mounting console with polymer shock absorbers (steel spring absorbers optionally)

HeNe - 633 - 0,7 - F - 117A - 2 - A - 28 - 0 - 150 **Order Code**

### Supply voltage:

230V (standard) ..... 2  
110V ..... 1

### Fiber length [cm]


### Fiber / connector type:

20 = PM singlemode fiber cable, FC-PC connector (0°-polish)  
28 = PM singlemode fiber cable, FC-APC connector (8°-polish), standard


### Fiber Connectors Type FC-APC

In order to avoid back reflection directly into the laser source, the fiber connector of type FC-APC has an 8° angled polish of the fiber ferrule. For fiber cables with this type of connectors Schäfter+Kirchhoff provides fiber collimators and laser beam couplers with inclined coupling axis.

### Dimensional drawings laser head



## Accessories Fiber Collimators 60FC-...



Fiber collimators made by Schäfter+ Kirchhoff transform radiation emitted divergently at the end of a singlemode fiber to a collimated beam. Standard collimators have a FC input connector, with inclined or coaxial coupling axis for FC-APC (angled polish) or FC-PC (plane polish) fiber connectors, respectively. Fiber collimators with ST, DIN AVIO or F-SMA are also available. On the front, the collimators have a receptacle to connect optics attachments such as micro focus optics for generation of micro spots down to 0.6 µm, or polarizers.

- Choice of 20 different collimating lenses, focal length ranging from 2.7 to 150 mm, for collimated beam diameters from 0.5 to 28 mm
- Various AR coatings in the range 370 - 2300 nm, with a bandwidth of a several 100 nm each.
- Lens focussing with indirect clamping, change of focus with attachments adapted


## Components for do-it-yourself laser coupling



A fiber coupling increases the possible field of application for existing HeNe lasers. The figure shows a 10 year old laboratory laser which was coupled to a polarization maintaining single-mode fiber with components from Schäfter + Kirchhoff. The necessary components, like laser beam coupler with

coupling optics adapted to the diameter of the incoming beam, mechanics adapters, and polarization maintaining singlemode fibers are provided separately. Additionally, Schäfter+Kirchhoff offers the fiber coupling of customer owned lasers as a service.


### Components from Schäfter+Kirchhoff



**Laser Beam Coupler 60SMS-...**  
Adjustable and focusable for singlemode fiber cable with FC connector  
Inclined or paraxial fiber coupling axis

The inclined fiber coupling axis is designed for fiber connectors of type FC-APC with 8° angled polish, and reducing back reflexion into the laser sources.

Download:  
[www.SuKHamburg.de/dl/lbc\\_e.pdf](http://www.SuKHamburg.de/dl/lbc_e.pdf)




#### Adapter Flange

For components with Ø 8 mm and Ø 19.5 mm system mount to adapt to Ø 19.5 mm mount.

#### Adapter Flange with Attenuator

Aperture Ø 3 mm with Ø 19.5 mm syst. mount




#### Polarization maintaining and Singlemode Fiber Cable PMC-.../ SMC-...

Singlemode and polarization maintaining, MFD 3 - 10 µm  
 $\lambda = 400 - 1800 \text{ nm}$

**Fiber connectors:**  
FC-APC: 8° angled polish of the fiber ferrule, suppression of back reflexion into the laser source.  
FC-PC: 0° polish

Download:  
[www.SuKHamburg.de/dl/fiber\\_e.pdf](http://www.SuKHamburg.de/dl/fiber_e.pdf)




#### Fiber Collimator 60FC... Focusable

Inclined or paraxial fiber coupling axis. Beam diameter and beam divergence are determined by the focal length  $f'$  of the collimating lens.


**Microfocus Optics, Series 5M-... and 13M-...**  
These lens attachments focus the collimated laser beam to a diffraction limited micro spot ( $> 0.001 \text{ mm}$ ).

Download:  
[www.SuKHamburg.de/dl/fc\\_e.pdf](http://www.SuKHamburg.de/dl/fc_e.pdf)


### Mounting consoles and accessories




**Mounting Console MC-MG-44,5-R** for HeNe lasers with diameter 44.5 mm/1.75" (set of two). Elastomer shock mounts are used for damping of shock and vibrations. Thermal deformation of the laser housing, which would be present in case of a rigid mounting of the laser, is prevented.  
The adapter 60A19,5-F or 60A19,5-F-AT is attached to the front plate of the laser.



**Mounting Console MC-MG-44,5-F-R** for HeNe lasers with diameter 44.5 mm/1.75" (set of two). Elastomer shock mounts for damping of shock and vibrations, and for avoidance of thermal deformations. For optimum stability, the console MC-MG-44,5-F-R holds both laser and adapter.



**Mounting Console MC-JDS-31,8-F-R** for HeNe lasers with diameter 31.8 mm/1.25" (set of two). Elastomer shock mounts for damping of shock and vibrations, and for avoidance of thermal deformations. For optimum stability, the console MC-MG-44,5-F-R holds both laser and adapter.



Socked head screws DIN912 4-40 UNC x 3/8", set of 4 pcs.

**Order Code**  
48-4-40-3/8-912-4


fitting tool: hex key  
**Order Code**  
50HD-3/32




## Dimensional Drawings for Power Supplies

HeNe lasers from Schäfter+Kirchhoff are supplied along with power supplies.  
Desktop power supplies are available for 230 V and 110 V line voltage.  
For some HeNe lasers there are OEM power supplies with input voltage 12 V DC or 230/100 V AC.


### Power Supply Style A



### Power Supply Style B




### Power Supply Style C




### OEM Power Supplies

By now, OEM power supplies with input voltage 12 V DC or 230/100 V AC are available for HeNe laser HeNe-633-0,6-P-1107P (Table 1 line 1). OEM power supplies for other HeNe lasers on request.

#### Power supply 115/230 V AC



#### Power supply 12 V DC



### Notes

## Safety at Work



### Laser safety and laser adjustment goggles

The use of laser safety goggles is recommended for work with lasers from laser protection class 3R on, e.g. all visible lasers from Schäfter+Kirchhoff with an output power up to 5 mW.

Laser safety goggles are mandatory from protection class 3B on, e.g. for all invisible infrared lasers and all visible lasers with more than 5 mW output power from Schäfter+Kirchhoff.

The correct handling and use of the laser safety goggles protects you and your colleagues against injuries of the eyes due to hazardous laser radiation.

A selection of CE and GS certified laser safety goggles is

available (manufactured by LaserVision, [www.lvg.com](http://www.lvg.com)) for the lasers provided by Schäfter+Kirchhoff.

Two different protective functions of laser safety goggles have to be distinguished, full protection goggles and alignment goggles.

Depending on filter type two different frames are supplied.

Laser safety goggles with glass filters have a frame for which an adapter frame (Order Code **RX7**) for individual glasses is available.

Laser safety goggles with plastic filter are more lightweight and can be worn together with glasses.

### Laser Safety Goggles - Function and Characteristics

**Protective function** Two different protective functions of laser safety goggles have to be distinguished, full protection goggles and alignment goggles.

**Full protection goggles** according to European standard EN 207 provide personal protection against laser radiation. The laser radiation normally is blocked and no longer visible. **Protection levels** (protection level L..) differ in the maximum spectral transmission of the filter glasses. In addition, the standard EN 207 specifies a maximum power density (power divided by area, in W/m²) which is allowed to reach the filter glass.

**Alignment protection goggles** according to European standard EN 208 attenuate visible laser radiation (wavelength 400 - 700 nm) down to the power of a laser of class 2 (EN 60825-1). The laser radiation remains visible, so that alignment protection glasses can be used for adjustment tasks. The **protection levels** (protection level R..) differ in the maximum power (in Watt) of a collimated laser beam which is allowed to reach the goggles.

**Maximum power** (EN 208): maximum power of a laser beam in the specified wavelength range which is sufficiently attenuated by the alignment protection goggles (acc. to EN 208).

**Maximum transmission** (EN 207): maximum trans-

mission (minimum attenuation) in the specified wavelength range (acc. to EN 208).

**Maximum power density** (EN 207): maximum power density the filter glasses can withstand over a longer period (acc. to EN 207)

**VLT**: (visible light transmission): in addition to the specified wavelengths, laser protection goggles also attenuate the ambient light. The VLT is the percentage of the daylight transmitted.

**OD (optical density)**: logarithmic scale for the attenuation of a laser beam at the specified wavelength. The OD at wavelength  $\lambda$  is defined as

$$OD(\lambda) = -\log_{10} \tau(\lambda)$$

### Goggles Types

#### Full and Alignment Protection Goggles acc. DIN EN 208/ DIN EN 207

Order Code **620.P0004.00**



VLT = 40%




### Usable Range

| Protection | Wave-length [nm] | Protection Level | max. Transmission (EN 207) | max. Power Density (EN 207) | max. Power (EN 208) |
|------------|------------------|------------------|----------------------------|-----------------------------|---------------------|
| Alignment  | 620- 644         | R2               | -                          | -                           | 100 mW              |
| Full       | 665- 685         | L4               | 10 <sup>-4</sup>           | 10 <sup>6</sup> W/m²        | -                   |


Alignment protection goggles for laser with wavelength 620-644 nm. For 665-685 nm full protection goggles.

### Optical Density OD



**Please note:** Shown are the typical density curves of the respective filters. These are no guaranteed values, only the protection levels (R.. or L..) are guaranteed by Schäfter+Kirchhoff.

## IR Detection Card





IR Detection Card SK100R Left: active sensor area Right: plot of spectral sensitivity

Invisible infrared radiation in the range 1100 - 1800 nm is converted to visible light by luminescence. The luminescence in the active detector layer indicates the position and approximate size of the beam.

The detector cards should not be used for measurement tasks like the quantitative determination of the beam diameter. With, e.g., high power lasers, also the border area of the beam with less than <1 % of the intensity maximum will be visible, additional widened by scattering processes in the detector layer. On the other hand, of a low power beam only the central area will be shown.

Detector cards of series SK ...R diffusely "reflect" the luminescent radiation at the active layer. These cards are impermeable for the beam. Detector cards of series SK...T are transparent. The beam is transmitted through the layer, the diffusely luminescent radiation is also visible from the back side of the card.



| Order Code         | Reflection SK 130R | Transmission SK 130T |
|--------------------|--------------------|----------------------|
| max. sensitivity   | 1300 nm            | 1300 nm              |
| Spectral range     | 1100 -1800 nm      | 1100 -1800 nm        |
| active sensor area | 20 x 20 mm²        | 20 x 20 mm²          |

## Laser Safety



**LASER LIGHT**  
**AVOID EXPOSURE TO BEAM**  
**CLASS III B LASER PRODUCT (CDRH)**  
**CLASS 3B LASER PRODUCT (DIN EN 60825-1-2001-11)**



## OPTICS, METROLOGY AND PHOTONICS BY

Schäfter + Kirchhoff GmbH

40 years of experience and a lot of production know-how are the base of many advanced optical products and optoelectronic systems designed by Schäfter + Kirchhoff

Laser Diode Systems in space

Optical Design Service

Customized Line Scan Camera

3D Object Profiling

Laser and Optics for AFM

Photolytic CNS stimulation

